

## 2-Color Wash Bottles—Safety-Labeled, Wide-Mouth Informative Labeling Available in 16 or 32 oz



### LAB-11646 Series

**Product Discontinued**

OMEGA's Safety Labeled wash bottles include a 53 mm (2.08") wide polypropylene closure with leakproof design for dispensing your most common solvents. Easy-to-read, LDPE solvent-resistant, printed bottles include: Department of Transportation (DOT) and National Fire Protection Association (NFPA) codes and symbols, Chemical Abstract Service (CAS) number, protective equipment symbols and target organ information, making your workplace a safer place. 6 per bag, 2 bags per case.

| Model No.      | Label                        | Cap Color | Capacity ml (oz) |
|----------------|------------------------------|-----------|------------------|
| LAB-116462232  | Acetone                      | Red       | 1000 (32)        |
| LAB-116466032  | Distilled Water              | Blue      | 1000 (32)        |
| LAB-116463739  | Ethanol                      | Natural   | 1000 (32)        |
| LAB-116462432  | Isopropanol                  | Yellow    | 1000 (32)        |
| LAB-116463832* | LYOB (Label Your Own Bottle) | Natural   | 1000 (32)        |
| LAB-116461132  | Methyl Ethyl Ketone (MEK)    | Green     | 1000 (32)        |

| Model No.      | Label                        | Cap Color | Capacity ml (oz) |
|----------------|------------------------------|-----------|------------------|
| LAB-116460050  | Assortment †                 | —         | 500 (16)         |
| LAB-116460622  | Acetone                      | Red       | 500 (16)         |
| LAB-116466250  | Dichloromethane              | Yellow    | 500 (16)         |
| LAB-116460631  | Deionized Water              | Blue      | 500 (16)         |
| LAB-116460620  | Distilled Water              | Blue      | 500 (16)         |
| LAB-116460639  | Ethanol                      | Natural   | 500 (16)         |
| LAB-116460640  | 70% Ethanol                  | Green     | 500 (16)         |
| LAB-116460612  | Ethyl Acetate                | Green     | 500 (16)         |
| LAB-116460624  | Isopropanol                  | Yellow    | 500 (16)         |
| LAB-116460638* | LYOB (Label Your Own Bottle) | Natural   | 500 (16)         |
| LAB-116460459  | Machine Oil (No Diamond)     | Natural   | 500 (16)         |
| LAB-116460623  | Methanol                     | Green     | 500 (16)         |
| LAB-116460611  | Methyl Ethyl Ketone (MEK)    | Green     | 500 (16)         |
| LAB-116460617  | Saline Solution (No Diamond) | Natural   | 500 (16)         |
| LAB-116460614  | Soap (No Diamond)            | Blue      | 500 (16)         |
| LAB-116460627  | Sodium Hypochlorite (Bleach) | Yellow    | 500 (16)         |
| LAB-116460628  | Toluene                      | Red       | 500 (16)         |
| LAB-116460621  | Water                        | Blue      | 500 (16)         |

\*Note: Use only with solvents compatible with LDPE.

† Assortment includes one each: Acetone, Ethanol, Isopropanol, Methanol, Toluene and Water.

Note: Sodium Hypochlorite bottle is white LDPE for light-protection between 230 and 450 nm. Toluene is red LDPE bottle.