




## Biohazard Disposal–Clear Bags Biohazard Label and Indicator Patch Option Available

### LAB-1316-Bag Series

#### Product Discontinued

Great economical solution for safe disposal of used pipettes, Petri dishes, multi-well cell culture plates, cell culture flasks, culture plates, membrane filters and other biohazard wastes. These polypropylene transparent biohazard bags are 0.03 mm (1.5 mil) thick.

Pass ASTM 1709 Dart® impact and ASTM 1922 tear resistance tests. High temperature autoclavable at a maximum of 135°C (275°F).



Available with or without Biohazard warning label.



Both models shown smaller than actual size.

#### A. Without Biohazard Warning Label or Sterilization Indicator Patch

Plain, non-printed bags can also be used as conventional trash can liners. 100 bags per box.

| To Order      |                   |                       |                 |
|---------------|-------------------|-----------------------|-----------------|
| Model No.     | Thickness         | Size                  | Gallon Capacity |
| LAB-131600005 | 0.03 mm (1.5 mil) | 30 x 61 cm (12 x 24") | 1 to 3          |
| LAB-131610005 | 0.03 mm (1.5 mil) | 61 x 76 cm (24 x 30") | 10 to 12        |
| LAB-131620005 | 0.03 mm (1.5 mil) | 61 x 91 cm (24 x 36") | 15 to 20        |

#### B. With Biohazard Warning Label and Sterilization Indicator Patch

Printed with the standard biohazard warning symbol and precautionary procedures in four languages: English, French, German, and Spanish. Sterilization patch darkens when subjected to steam sterilization. 100 bags per box.

| Model No.     | Thickness | Size                  | Gallon Capacity |
|---------------|-----------|-----------------------|-----------------|
| LAB-131600009 | 1.5 mil   | 30 x 61 cm (12 x 24") | 1 to 3          |
| LAB-131610009 | 1.5 mil   | 61 x 76 cm (24 x 30") | 10 to 12        |
| LAB-131620009 | 1.5 mil   | 61 x 91 cm (24 x 36") | 15 to 20        |