

## Benchtop Biohazard Bags Easy Open, Dispose Used Tips and Tubs

### LAB-13166 Series

Sold per box, 100 or 1000

**Product Discontinued**

Bright red, autoclavable bags feature one long flap for easy opening. Their size allows the bag to be inverted over your hand to safely pick up waste. For a space-saving disposal unit in fume hoods or on bench tops, use them with Poxgrid® desk model bag holder LAB-131931000. Printed with the biohazard symbol and sterilization indicator patch, they are made of high molecular weight, high-density (HMHD) polyethylene and pass ASTM 1922 tear resistance and ASTM 1709 Dart® impact tests. Steam autoclavable at 121°C (250°F). Sold per box, 100 or 1000.



Shown smaller than actual size.

### Quick Method for Waste Pick-up

*Insert gloved hand into bag.*

*Pick up waste materials with hand inside bag.*

*Peel bag over the hand, turning it inside out while maintaining grip on waste material.*

*Waste is now inside bag and ready for autoclaving.*

To Order				
Model No.	Size	Quantity	Description	Thickness
LAB-131660000	0.43 gallon	100	Benchtop biohazard bag, 22 x 28 cm H (8.5 x 11")	0.018 mm (0.72 mil)
LAB-131660001	0.43 gallon	1000	Benchtop biohazard bag, 22 x 28 cm H (8.5 x 11")	0.018 mm (0.72 mil)