

Multipurpose Trays

LAB-16200 Series

Product Discontinued



LAB-162000010
multipurpose tray
shown smaller
than actual size.



LAB-162000011
multipurpose tray with
handles shown smaller
than actual size.

These large capacity, rugged polypropylene trays have multiple uses. The substantial depth and 8.5 liter (9 quart) volume of each can be used to hold a large quantity of samples in a water bath or on ice. The trays are steam autoclavable at 121°C (250°F) and ideal for holding instruments while autoclaving. Use in a sink when cleaning glassware to reduce chance of damage or for transporting labware. Inside corners are rounded for easy cleaning. The top rim adds rigidity and provides handholds for transport. Each size nests with identical trays to minimize storage space requirements.

To Order	
Model No.	Inside Dimensions cm (inch)
LAB-162000010	30.5 L x 24.8 W x 11.4 H (12 x 9¾ x 4½)
LAB-162000011	32.4 L x 26.7 W x 10.8 H (12¾ x 10½ x 4¼)