

Clips for Beakers and Pitchers




LAB-182350000 Series

Product Discontinued

- Sturdy Holder Clips Onto The Side of an Open Container with a Wall up to 6.35 mm (¼") Thick
- Holder Has Room for Sample Identification Marks

These easy-to-use, polypropylene clips affix to the rim of a beaker or pitcher to hold tubing, dialysis bags or thermometers.

All models shown larger than actual size.



| To Order | | | |
|---------------|-------------------------|----------------------------------|------------------|
| Model No. | Description | Holds | Quantity per Bag |
| LAB-182350000 | Tubing holder "A" | Up to 9.5 mm (3/8") O.D. | 6 |
| LAB-182370000 | Dialysis bag holder "B" | On rim up to ¼" thick | 6 |
| LAB-182410000 | Thermometer holder "C" | 3.2 to 11 mm dia. (1/8 to 7/16") | 3 |