

Large Weighted Safety Shields
Highest Impact Resistance of All Transparent Shield Materials

LAB-24964 Series

Laboratory workers are protected by a 4.7 mm ($\frac{3}{16}$ ") thick clear polycarbonate barrier. Stable three-sided construction provides front and side protection against flying debris. The heavy steel perimeter base is coated with blue epoxy. Plated steel bolts, nuts and washers are used for assembly. Outside base width is 50 cm (19 $\frac{1}{2}$ ") and 17 cm (6 $\frac{1}{2}$ ") front to back.

Product Discontinued



LAB-249640000 shown smaller than actual size.

To Order		
Model No.	Height cm (inch)	Weight kg (lb)
LAB-249620000	76 (30)	8.4 (18.5)
LAB-249640000	91 (36)	8.9 (19.6)