

Beta Shields
Free Standing Workshields in 3 Sizes

LAB-24976 Series

Product Discontinued

Beta shields reduce the exposure from radioactive isotopes and offer excellent clarity. The shields are made from 9.5 mm (3/8") clear acrylic with non-skid feet. All three sizes feature an angled shield and flat base with large work surface.



To Order		
Model No.	Description	Dimensions cm (inch)
LAB-249760001	Acrylic free standing beta shield	20.3 L x 30.5 W x 37.5 H (8 x 12 x 14 ³ / ₄)
LAB-249760002	Acrylic free standing beta shield	30.5 L x 35.6 W x 52.7 H (12 x 14 x 20 ³ / ₄)
LAB-249760003	Acrylic free standing beta shield	31.4 L x 40.6 W x 52.1 H (12 ³ / ₈ x 16 x 20 ¹ / ₂)