

## Graduated Cylinders with Handles

**Product Discontinued**



### LAB-28461 Series

- Grip Handle for Steady Pouring

This one-piece polypropylene cylinder with handle provides a secure, solid grip for safer pouring. The “total contained, total delivered” double scale is easy-to-read and accurate to  $\pm 5\%$ . A hexagonal base prevents tipping and rolling. Handle opening is 8.9 x 2.9 cm (3 x 1 $\frac{1}{8}$ ”). Chemically resistant and autoclavable at 121°C (250°F).

To Order					
Model No.	Capacity	Graduation	Height	I.D.	O.D.
LAB-284611000	1000 ml	Every 10 ml	439 mm (17.3")	58 mm (2.3")	64 mm (2.5")
LAB-284612000	2000 ml	Every 20 ml	490 mm (19.3")	84 mm (3.3")	89 mm (3.5")