



Polypropylene Scoops Six Sizes

Product Discontinued

LAB-367-SCOOP Series

These scoops are rigid, strong and smooth. A flat bottom keeps them upright when filled or empty, allowing them to be used as weighing boats or to be placed on a flat surface. The 60 and 12, 5 ml scoops have a longer handle that increases the distance between the hand and sample. Steam autoclavable at 121°C (250°F), they comply with FDA regulations for food contact. Special prices on yearly contracts for large quantities and colors are available on request. 12 per bag except *LAB-36756000 (6 per bag). Bulk quantity pricing available.



Shown smaller than actual size.

To Order				
Model No.	Capacity ml (oz)	Bowl (L x W x H)	Handle	Package Quantity
LAB-367470000	22 (0.74)	6 x 2.9 x 2.5 cm (2 ³ / ₈ x 1 ¹ / ₈ x 1")	3.2 cm (1 ¹ / ₄ ")	12
LAB-367480000	60 (2)	7.6 x 3.8 x 2.5 cm (3 x 1 ¹ / ₂ x 1")	15 cm (6")	12
LAB-367490000	75 (2.5)	8.9 x 3.8 x 3.2 cm (3 ¹ / ₂ x 1 ¹ / ₂ x 1 ¹ / ₄ ")	4.8 cm (1 ⁷ / ₈ ")	12
LAB-367500000	125 (4)	9.2 x 4.8 x 3.5 cm (3 ⁵ / ₈ x 1 ⁷ / ₈ x 1 ³ / ₈ ")	15 cm (6")	12
LAB-367530000	250 (8)	12.7 x 5.1 x 5.1 cm (5 x 2 x 2")	6.4 cm (2 ¹ / ₂ ")	12
LAB-367560000*	1100 (37.2)	20.3 x 12.1 x 7.6 cm (8 x 4 ³ / ₄ x 3")	10.2 cm (4")	6

These sampling scoops are designed for retrieving samples from a larger source. The capacity/volumes are an estimated volume and are not meant for taking exact sample weights/volumes. Stated volumes represent typical heaping samples.