

Liquid Nitrogen Cooled Mortar
Grind Frozen Samples Without Liquid Nitrogen Hazards

LAB-37260 Series

Product Discontinued

A 304 stainless steel mortar is suspended by three brackets above a stainless steel bowl which holds liquid nitrogen. The bowl is nested in a heavily insulated blue polyethylene housing. The mortar has an insulated handle, 10.2 cm L (4"), that is removable for saving space. The housing measures 16.5 cm dia (6½") x 11.4 cm high (4½") and has three suction cup feet. The recommended size pestle is any size appropriate for a 275 ml mortar. Pestle not included.



To Order	
Model No.	Description
LAB-372600000	Stainless steel ladle (bowl and detachable handle) and stainless steel reservoir in a polyethylene housing
LAB-372600001	Spare stainless steel mortar