

**Pipette Aid**  
**Safe, Easy and Economical Pipetting**

**LAB-378860000**

**Product Discontinued**

This pipette aid is used by squeezing, then releasing the bulb to draw fluid into the pipette. The pipette aid is then removed, and the pipette top quickly covered with a gloved finger to control and release fluid. The bulb is Vikem® vinyl with a tapered polyethylene chuck for quick removal when pipette is filled. 6 per pack.

**Dimensions:** 5 dia x 10 cm H (1.9 x 3.9")



**To Order**

Model No.	Description
LAB-378860000	Pipette aid