

Pipette Pump Pipettors Four Sizes

LAB-3789 Series

Product Discontinued

Provides safe, accurate and trouble-free pipetting with simple, one-handed operation. Pipettes fit smoothly into its flexible tapered chuck. The thumbwheel on the side is rotated for precision filling or dispensing, and the plunger may be pressed for quick emptying. It is resistant to acids and alkalis. Disassemble all parts for cleaning and before autoclaving at 121°C (250°F) for 20 minutes except the chuck which should not be autoclaved. Sold individually.



All models shown smaller than actual size.

To Order		
Model No.	Pipette Up To	Color
LAB-378960000	0.2 ml	Yellow
LAB-378970000	2 ml	Blue
LAB-378980000	10 ml	Green
LAB-378990000	25 ml	Red

Replacement Silicone Rubber Chuck and Washers

Model No.	Pipette Chuck	Washers
LAB-378960000	LAB-937890060	LAB-937890061 sleeve
LAB-378970000	LAB-937890070	LAB-937890071, 3/Set
LAB-378980000	LAB-937890080	LAB-937890081, 3/Set
LAB-378990000	LAB-937890090	LAB-937890091, 3/Set