



Economy Pipette Pump III ***Accuracy at an Affordable Price!***

LAB-37902 Series

Product Discontinued

- Easy One-Handed Operation
- 2, 10, and 25 ml Pipette Volumes
- 3 Color Codes for Easy Identification

Ergonomic design enhances comfort, utility and dispensing precision while helping to reduce the risk of injury from repetitive and awkward pipetting motions. This lightweight, affordable device provides accurate, trouble-free pipetting with simple, one-handed operation. The optimally located thumb wheel rotates easily for precision aspirating or dispensing, down to a volume of one drop at a time. The thumb wheel also can be used to "blow out" the contents of the pipette. Available in three different pipette volumes, 2, 10 and 25 ml, each color-coded for easy identification. Resistant to acids, alkalis and most solvents, they can be disassembled for cleaning. Sold individually.

Overall Length: 17.1 cm (6¾")



To Order

Model No.	Pipette Up To	Color
LAB-379020002	2 ml	Blue
LAB-379020010	10 ml	Green
LAB-379020025	25 ml	Red