



**Fast Release Pipette Pump**  
**Features Rapid Release Valve for Quick Dispensing**

**Product Discontinued**

## LAB-37911 Series

Glass or plastic pipettes are easily inserted into a unique collar with stabilizing fingers for a secure hold. Rotating the knurled thumb wheel draws liquid to the desired volume. To empty rapidly, apply light pressure to the fast release lever. For gradual dispensing, rotate the thumb wheel. Sizes are color coded. They can be disassembled for cleaning. Sold individually.



All models shown  
smaller than actual size.

To Order		
Model No.	Pipette Up To	Color
LAB-379111002	2 ml	Blue
LAB-379111010	10 ml	Green
LAB-379111025	25 ml	Red