

Ampule Opener

LAB-389520000

Product Discontinued

- Includes 4 Inserts to Allow Opening of 1, 2, 5 and 10 ml Ampules
- Open Ampules Easily and Safely
- Transparency Allows Unobstructed Viewing of Ampule Contents at all Times
- Copolyester Plastic Construction has Broad Chemical Resistance

Box style design allows one handed opening of up to 3 ampules at a time by simply pressing down on the lid. Ampule necks snap clean and are contained within the device while fingers stay away from sharp edges. After removing ampule contents, turn the device upside down over a trash receptacle to safely discard ampules and tops.

Dimensions: 9.5 L x 7.6 W x 8.3 cm H (3.75 x 3 x 3.25")



LAB-389520000 shown smaller than actual size.

To Order	
Model No.	Description
LAB-389520000	Ampule opener