

Vacuum Chamber and Plate *For Vacuum Drying, Degassing and Teaching Demonstrations*

LAB-420430000

Product Discontinued

- Interior Volume: 0.21 Cu. Ft.
- Can Withstand a Full Vacuum at Room Temperature
- Clear Polycarbonate Top
- ABS Plastic Plate and O-Ring Seals the Top Against the Chamber's Flange
- Stopcock with 2 mm (0.08") PTFE Plug
Uses 6 mm (¼") Tubing

LAB-420430000 shown smaller than actual size.



To Order

Model No.	Description	Dimensions (Diameter x H)
LAB-420430000	Vacuum chamber and plate	22.8 x 16 cm (9 x 6¼")