



## **Small Molded Benchtop Fume Hood Saves Space and Money**

**LAB-500000002**

**Product Discontinued**

Designed for small facilities or installations in which a number of fume hoods are needed, this unit is excellent for use in schools, junior colleges and universities. Features include a one-piece shell molded of chemical resistant polyethylene with rounded inside corners for easy cleaning. The clear acrylic sash slides on a track featuring stop pins for open, closed or intermediate positions, and a bottom gap allows for continuous air flow. The stepped exhaust duct can be used with either a 12 or 15 cm (4 or 6") diameter duct.

**Overall Dimensions:** 69 L x 33 W x 56 cm H (27 x 13 x 22")  
**Front Opening:** 56 W x 30 cm H (22 x 12")  
**Exhaust System Requirements:** 200 cfm (free air) recommended  
**Shipping Weight:** 17 kg (38 lb)



### To Order

Model No.	Description	Dimensions (L x W x H)
LAB-500000002	Small molded benchtop fume hood	69 x 33 x 56 cm (27 x 13 x 22")