

Observation Fume Hood **Maximum Visibility and Protection During Lessons**

LAB-500002002

Lightweight and economical, this hood has clear polycarbonate windows which are stronger than acrylic, wipe clean with a damp cloth, and provide excellent visualization in a group setting. The front sash opens to three adjustable positions for air flow control. The hood's polyethylene surface is chemical resistant and easy to clean. Stepped exhaust duct will accept either a 12 or 15 cm (4 or 6") diameter duct.

Product Discontinued

Overall Dimensions: 69 L x 33 W x 56 cm H (27 x 13 x 22")

Back Stationary Window: 50.8 W x 27.9 cm H (20 x 11")

Exhaust System Requirements: 200 cfm (free air) recommended

Shipping Weight: 18.1 kg (40 lb)



To Order

Model No.	Description	Dimensions (L x W x H)
LAB-500002002	Observation fume hood	69 x 33 x 56 cm (27 x 13 x 22")