



## Portable pH Meters



**PHH-65D**

**1 YEAR  
ELECTRONICS  
WARRANTY**

OMEGA® portable pH meters are designed to be reliable, low-cost instrument for pH, specific ion, and redox potential measurements directly at the sampling site. The PHH-65A and PHH-65D are equipped with 8 AA batteries and are housed in rugged, high-impact cases that contain and protects the meter and associated accessories. For ease of use, simple operating instructions are contained in the inside cover for handy reference. These units come complete with a gel-filled combination pH electrode, buffer solutions (pH 4, 7, 10) and operator's manual.

### PHH-65D Digital pH Meter

The PHH-65D meter features a 3½-digit, 12.7 mm (½") high LCD display that is easy to read even in the brightest sunlight. This instrument is capable of measuring pH values from 0 to 14 with an accuracy of 0.05 pH and 0 to 1,999 mV with an accuracy of ±1 mV.



PHH-65A complete kit as shown.

### PHH-65A Complete Kit

#### PHH-65A Analog pH Meter

The PHH-65A meter features a large 6-inch analog scale for reading both normal (pH 0 to 14, ±700 mV) and 10x expanded scale reading (±70 mV and any 1.4 pH full scale). Using the expanded scale, accuracies as great as 0.01 pH can be achieved for critical applications.

#### Specifications

- Range:** 0 to 14 pH, ±700 mV (±1999 mV with PHH-65D)
- Accuracy:** ±0.05 pH, PHH-65A; ±0.01 pH, PHH-65D and PHH-65A expanded scale; ±1 mV, both models
- Input Impedance:** >10<sup>10</sup> Ω
- Temperature Compensation:** 0 to 100°C (32 to 212°F), manual
- Connector:** BNC
- Power:** 8 "AA" batteries (included)
- Dimensions:** 114 H x 311 W x 203 mm D (4½ x 12¼ x 8")
- Weight:** 1.8 kg (4 lb)
- Optional Accessories:** Recommended replacement electrode, PHE-1304