

Compact Stirrer

PHS-21

- Small Yet Powerful, 127 x 127 mm (5 x 5") Size

The PHS-21 is a compact, efficient and economical general purpose stirrer. It has been developed to fit the needs of today's lab.

Small yet powerful, the 127 x 127 mm (5 x 5") size fits neatly into nearly any lab space and the Brushless 2 W motor will give years of trouble free service. The bright stainless steel top provides better stirring power and is easy to clean.



SPECIFICATIONS

Top Plate Material: Bright stainless steel
Stirrer Speed: 300 to 2800 RPM
Maximum Stirring Capacity: 2.8 liters
Motor: 2 W, brushless
Dimension: 128 W x 128 L x 55 D mm (5 x 5 x 2.1")
Plate Dimensions: 127 x 127 mm (5 x 5")
Net Weight: 0.8 kg (1.7 lb)
Power Supply: AC 90 to 230V, 50 to 60 Hz

To Order

Model No.	Description
PHS-21	Stainless steel stirrer with universal power output analog